

Date:

Name veterinary practice:

Practice address:

Postcode and city:

Country:

Artuvetrin® Skin Test order form

Select your choice of allergens

This is a selection of the most frequent occurring allergens. Please consult our product range for additional allergens available or contact our sales department: +31 (0)320 783 100.

Basic Skin Test Set

☐ Artuvetrin® Skin Test Set

- Physiological phosphate buffer (Negative control)
- Histamine phosphate (Positive control)
- Tree pollen mixture 1*
- Tree pollen mixture 2*
- Weed pollen mixture*
- Grass pollen mixture*
- Fungi mixture 1*
- Acarus siro (Grain mite)
- Derm. farinae (Farinae mite)
- Derm. pteronyssinus (House dust mite)
- Lepidoglyphus destructor (Hay mite)
- Tyrophagus putrescentiae (Copra mite)
- Felis domesticus (Cat-epithelia)

Individual Skin Tests

- ☐ Physiological phosphate buffer (Negative control)
- ☐ Histamine phosphate (positive control)

Grass pollen

- ☐ Grass pollen mixture*
- ☐ Anthoxanthum odoratum (Sweet vernal grass)*
- ☐ Festuca pratensis (Meadow fescue)
- ☐ Holcus lanatus (Velvet grass)*
- ☐ Lolium multiflorum (Italian rye grass)
- ☐ Lolium perenne (Perennial rye grass)
- ☐ Phleum pratense (Timothy)*
- ☐ Poa pratensis (Kentucky blue grass)

Pollen of crops

- ☐ Avena sativa (Cultivated oat)
- ☐ Brassica napus (Rape)
- ☐ Secale cereale (Common rye)
- ☐ Triticum aestivum (Wheat)
- ☐ Zea mays (Corn)

Pollen from weeds and flowers

- ☐ Weed pollen mixture*
- ☐ Weed pollen mixture 1*
- ☐ Weed pollen mixture 2*
- ☐ Artemisia vulgaris (Common mugwort)*
- ☐ Chenopodium album (Lamb's quarters)*
- ☐ Leucanthemum vulgare (Daisy)
- ☐ Plantago lanceolata (English plantain)*
- ☐ Rumex acetosella (Red or sheep sorrel)*
- ☐ Rumex crispus (Yellow dock)

Pollen from trees and hedges

- ☐ Tree pollen mixture 1*
- ☐ Tree pollen mixture 2*
- ☐ Olea Europea (Olive)
- ☐ Platanus occidentalis (Eastern sycamore)

Mites

- ☐ Acarus siro (Grain mite)
- ☐ Euroglyphus maynei
- ☐ Lepidoglyphus destructor (Hay mite)
- ☐ Derm. farinae (Farinae mite)
- ☐ Derm. pteronyssinus (House dust mite)
- ☐ Tyrophagus putrescentiae (Copra mite)

Fungi

- ☐ Fungi mixture 1*
- ☐ Aspergillus mix*
- ☐ Alternaria alternata*
- ☐ Aspergillus fumigatus*
- ☐ Cladosporium herbarum*

Insects

- ☐ Aedes communis (Mosquito)
- ☐ Culex (Mosquito)
- ☐ Culicoides (Biting midges)
- ☐ Musca domestica (House fly)
- ☐ Periplaneta Americana (Cockroach)
- ☐ Tabanus (Horse fly)

Epithelia

- ☐ Epithelia mixture 1*
- ☐ Epithelia mixture 2*
- ☐ Equus caballus (Horse)
- ☐ Felis domesticus (Cat)*
- ☐ Melopsittacus undulatus (Parakeet)
- ☐ Ovis aries (Sheep)

Yeast

- ☐ Malassezia

Other allergens

Others (please indicate number)

- ☐ Materials for Next+® Serum Test
- ☐ Order form Artuvetrin® Skin Test
- ☐ Order form Artuvetrin® Therapy
- ☐ Flow chart Allergic Dermatitis
- ☐ Product range Artuvetrin®
- ☐ Client brochure 'Allergy in dogs'
- ☐ Wheal measuring ruler
- ☐ Skin marker black
- ☐ Skin marker blue

*For more information about the allergen mixtures p.t.o.

Explanation allergen mixtures

Selection of the most frequent occurring allergens

Grass pollen mixture

- Bermuda grass
- Orchard grass
- Sweet vernal grass
- Timothy grass
- Velvet grass

Tree pollen mixture I

- Birch
- Alder
- Hazel

Tree pollen mixture 2

- English oak
- European beech
- American elm

Weed pollen mixture

- Common mugwort
- Stinging nettle
- Common dandelion
- English plantain

Weed pollen mixture I

- Common mugwort
- Stinging nettle

Weed pollen mixture 2

- Goldenrod
- Lamb's quarters
- Common dandelion
- Red or sheep sorrel
- English plantain

Fungi mixture 1

- Alternaria alternata
- Aspergillus fumigatus
- Cladosporium herbarum

Epithelia mixture 1

- Guinea pig
- Dog
- Cat
- Hamster
- Rabbit

Epithelia mixture 2

- Duck
- Goose
- Chicken

Aspergillus mixture

- Aspergillus flavus
- Aspergillus fumigatus
- Aspergillus nidulans
- Aspergillus niger

Note: This order form contains the most frequent occurring allergens. It is also possible to order other allergens. Information about all the different allergens supplied by Nextmune can be found in our product range brochure. This brochure as well as the order form Artuvetrin® Therapy can be ordered with this form, or by phone +31 (0)320 783 100 or email info.eu@nextmune.com.

Additional allergens can be filled out under 'Other allergens' (individual tests). Be aware that in order to do a Skin Test you always need a positive and negative control and that these tests need to be ordered separately with individual tests.

The Next+® Serum Test is a good and reliable alternative for the Artuvetrin® Skin Test. Information about this serological test can also be ordered with this form.